Drew Scientific Co. Ltd.
Glycated Haemoglobin Two Level Control Kit

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For in vitro diagnostic use only. These controls are intended for use in a quality control procedure for determining imprecision and accuracy of glycated haemoglobin measurements.

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The routine use of Control Materials is recommended for the long term monitoring of imprecision of analytical methods and is a requirement of Good Laboratory Practice.

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Contents, Contenu du kit, Inhalt, Contenuto, Conteúdo, Contenido

The kit contains three vials of each of two levels. Each vial contains a preparation of freeze dried whole blood.

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Storage, Conservation, Lagerung, Conservazione, Armazenagem, Almacenamiento

Store unopened at 2 – 8 °C. The expiration date is given at the top of this page. After reconstitution, the material is stable for 7 days at 2 – 8 °C.

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Manufactured by Drew Scientific
Peter Green Way, Barrow-in-Furness,
Cumbria LA14 2PE
Tel: +44(0)1229 432089 / Fax: +44(0)1229 432096

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Carefully remove stopper. Using a volumetric pipette, measure 0.3mL of deionised water and add to the vial to reconstitute the control. Replace stopper and gently mix. Allow to stand for 5 minutes. Visually inspect that all the material has dissolved before use. Treat the reconstituted control as whole blood for the purpose of preparing samples for your glycated haemoglobin analysis method.

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Do not use time-expired material. Do not use incompletely dissolved material.

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All the individual blood donor lots used to prepare this material were tested by FDA approved methods and found to be non-reactive for hepatitis B surface antigen (HbsAG), hepatitis C (HCV), HIV-1 and HIV-2. No test methods can offer complete assurance that these materials are free from infective agents. Handle, use and dispose of these control materials as if they were human specimens and potentially infective.

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Assigned values

Methods	LEVEL 1 CONTROL	
	HbA1c Value	Guidance Range
Drew Scientific DS5 <i>NGSP Aligned</i>	5.4 %	4.3 % to 6.5 %
Drew Scientific DS5 <i>IFCC Aligned</i>	3.6 % <i>(equivalent to 36 mmol/mol Hb)</i>	2.5 % to 4.7 %
Drew Scientific G15 <i>NGSP Aligned</i>	5.7 %	4.9 % to 6.5 %
Drew Scientific G15 <i>IFCC Aligned</i>	3.9 % <i>(equivalent to 39 mmol/mol Hb)</i>	3.1 % to 4.7 %

LEVEL 2 CONTROL	
HbA1c Value	Guidance Range
10.1 %	8.4 % to 11.8 %
8.6 % <i>(equivalent to 86 mmol/mol Hb)</i>	6.8 % to 10.4 %
11.0 %	9.9 % to 12.1 %
9.6 % <i>(equivalent to 96 mmol/mol Hb)</i>	8.4 % to 10.8 %