



Product No. 610020
Lot No. CDC0917
Expiration : 2017-10-05

HEMATOLOGY CALIBRATOR HEMAVET SERIES

For use with HEMAVET[®] HEMATOLOGY INSTRUMENTS

Parameter (units)	Target	Range
WBC (k/ μ L)	8.15	7.64 – 8.66
RBC (M/ μ L)	4.47	4.19 – 4.75
Hb (g/dL)	13.5	12.9 – 14.1
MCV (fL)	78.5	77.0 – 80.1
PLT (k/ μ L)	270	253 – 287

IMPORTANT

Follow mixing instructions from reverse side and
Wipe threads on bottle and cap before closing

**FOR IN VITRO DIAGNOSTIC USE
PROTECT FROM FREEZING
STORE AT 2°C - 8°C (35°F - 46°F)**

POTENTIAL BIOHAZARDOUS MATERIAL

DREW SCIENTIFIC, INC. 14100 NW 57th Court, Miami Lakes, FL 33014 • USA
Phone: (305) 418-2320 Fax: (305) 418-2321
Website: drew-scientific.com