



DREW SCIENTIFIC, INC.

MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: CAT Product Number: 600064, 600067, 600069

Lot 2004 Expiration Date: 10/15/2017

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 550,1500,1700,3700, Mascot 850		Hemavet 950, 1700 Rev F and above	
	Target	Range	Target	Range
WBC (K/μL)	7.1	6.5 - 7.7	7.1	6.5 - 7.7
NE (K/μL)				
NE/MO/BA (K/μL)				
LY (K/μL)	1.8	1.3 - 2.3	1.8	1.3 - 2.3
MO (K/μL)				--Elev
EO (K/μL)				
BA (K/μL)				--Elev
NRBC (K/μL)				
NE (%)				
NE/MO/BA (%)				
LY (%)	25.4	20.4 - 30.4	25.4	20.4 - 30.4
MO (%)				--Elev
EO (%)				
BA (%)				--Elev
NRBC (%)				
RBC (M/μL)	7.3	6.98 - 7.62	7.30	6.98 - 7.62
HB (g/dL)	13.3	12.7 - 13.9	13.3	12.7 - 13.9
HCT (%)	39.0	34.5 - 43.8	39.0	34.5 - 43.8
MCV (fL)	53.4	49.4 - 57.4	53.4	49.4 - 57.4
MCH (pg)	18.2	16.7 - 19.9	18.2	16.7 - 19.9
MCHC (g/dL)	34.1	29.0 - 40.3	34.1	29.0 - 40.3
RDW (%)	21.9	19.9 - 23.9	21.9	19.9 - 23.9
PLT (K/μL)	182	137 - 227	182	137 - 227

DREW SCIENTIFIC, INC. 14100 NW 57th Court Miami Lakes, FL 33014 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2°C - 8°C (35°F - 46°F)