



# DREW SCIENTIFIC, INC.

## MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: CAT      Product Number: 600064, 600067, 600069

Lot 2014      Expiration Date: 1/15/2018

**For use with HEMAVET® and MASCOT™ Hematology Analyzers**

Parameter	Hemavet 800,850,850FS,950,950LV 550,1500,1700,3700, Mascot 850		Hemavet 950, 1700 Rev F and above	
	Target	Range	Target	Range
WBC (K/μL)	7.1	6.5 - 7.7	7.1	6.5 - 7.7
NE (K/μL)				
NE/MO/BA (K/μL)				
LY (K/μL)	1.6	1.2 - 2.2	1.6	1.2 - 2.2
MO (K/μL)				--Elev
EO (K/μL)				
BA (K/μL)				--Elev
NRBC (K/μL)				
NE (%)				
NE/MO/BA (%)				
LY (%)	23.0	18.0 - 28.0	23.0	18.0 - 28.0
MO (%)				--Elev
EO (%)				
BA (%)				--Elev
NRBC (%)				
RBC (M/μL)	7.7	7.36 - 8.00	7.68	7.36 - 8.00
HB (g/dL)	13.3	12.7 - 13.9	13.3	12.7 - 13.9
HCT (%)	38.3	33.8 - 43.1	38.3	33.8 - 43.1
MCV (fL)	49.9	45.9 - 53.9	49.9	45.9 - 53.9
MCH (pg)	17.4	15.9 - 18.9	17.4	15.9 - 18.9
MCHC (g/dL)	34.9	29.6 - 41.3	34.8	29.6 - 41.3
RDW (%)	26.4	24.4 - 28.4	26.4	24.4 - 28.4
PLT (K/μL)	179	134 - 224	179	134 - 224

DREW SCIENTIFIC, INC. 14100 NW 57th Court Miami Lakes, FL 33014 • USA

### IMPORTANT

Follow mixing instructions from reverse side:  
Wipe threads on bottle and cap before closing

**FOR INVITRO DIAGNOSTIC USE**

**PROTECT FROM FREEZING**

STORE AT 2°C - 8°C (35°F - 46°F)