



DREW SCIENTIFIC, INC. MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: CAT Product Number: 600064, 600067, 600069

Lot 2034 Expiration Date: 7/15/2018

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 550,1500,1700,3700, Mascot 850		Hemavet 950, 1700 Rev F and above	
	Target	Range	Target	Range
WBC (K/μL)	7.5	6.9 - 8.1	7.5	6.9 - 8.1
NE (K/μL)	0.1			
NE/MO/BA (K/μL)				
LY (K/μL)	2.6	2.1 - 3.3	2.6	2.1 - 3.3
MO (K/μL)				--Elev
EO (K/μL)				
BA (K/μL)				--Elev
NRBC (K/μL)				
NE (%)				
NE/MO/BA (%)				
LY (%)	35.1	30.1 - 40.1	35.1	30.1 - 40.1
MO (%)				--Elev
EO (%)				
BA (%)				--Elev
NRBC (%)				
RBC (M/μL)	7.1	6.82 - 7.46	7.14	6.82 - 7.46
HB (g/dL)	12.7	12.1 - 13.3	12.7	12.1 - 13.3
HCT (%)	39.9	35.4 - 44.7	39.9	35.4 - 44.7
MCV (fL)	55.9	51.9 - 59.9	55.9	51.9 - 59.9
MCH (pg)	17.8	16.2 - 19.5	17.8	16.2 - 19.5
MCHC (g/dL)	31.8	27.0 - 37.5	31.8	27.0 - 37.5
RDW (%)	22.5	20.5 - 24.5	22.5	20.5 - 24.5
PLT (K/μL)	188	143 - 233	188	143 - 233

DREW SCIENTIFIC, INC. 14100 NW 57th Court Miami Lakes, FL 33014 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2°C - 8°C (35°F - 46°F)