



DREW SCIENTIFIC, INC. MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: CAT Product Number: 600064, 600067, 600069

Lot 2044 Expiration Date: 10/15/2018

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 550,1500,1700,3700, Mascot 850		Hemavet 950, 1700 Rev F and above	
	Target	Range	Target	Range
WBC (K/μL)	7.5	6.9 - 8.1	7.5	6.9 - 8.1
NE (K/μL)	0.1			
NE/MO/BA (K/μL)				
LY (K/μL)	1.3	0.9 - 1.9	1.3	0.9 - 1.9
MO (K/μL)				--Elev
EO (K/μL)				
BA (K/μL)				--Elev
NRBC (K/μL)				
NE (%)				
NE/MO/BA (%)				
LY (%)	17.8	12.8 - 22.8	17.8	12.8 - 22.8
MO (%)				--Elev
EO (%)				
BA (%)				--Elev
NRBC (%)				
RBC (M/μL)	7.3	6.95 - 7.59	7.27	6.95 - 7.59
HB (g/dL)	12.9	12.3 - 13.5	12.9	12.3 - 13.5
HCT (%)	40.1	35.6 - 44.9	40.1	35.6 - 44.9
MCV (fL)	55.2	51.2 - 59.2	55.2	51.2 - 59.2
MCH (pg)	17.7	16.2 - 19.4	17.7	16.2 - 19.4
MCHC (g/dL)	32.1	27.3 - 37.8	32.1	27.3 - 37.8
RDW (%)	22.7	20.7 - 24.7	22.7	20.7 - 24.7
PLT (K/μL)	182	137 - 227	182	137 - 227

DREW SCIENTIFIC, INC. 14100 NW 57th Court Miami Lakes, FL 33014 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2°C - 8°C (35°F - 46°F)