



# DREW SCIENTIFIC, INC.

## MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: DOG      Product Number: 600066, 600067, 600068

Lot 2016      Expiration Date: 1/15/2018

**For use with HEMAVET® and MASCOT™ Hematology Analyzers**

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	11.3	10.7 - 11.9	11.3	10.7 - 11.9
NE (K/μL)	7.5	6.2 - 8.8		
NE/MO/BA (K/μL)			8.3	7.0 - 9.7
LY (K/μL)	1.7	1.1 - 2.4	1.7	1.1 - 2.4
MO (K/μL)	0.8	0.0 - 1.1		--Elev
EO (K/μL)	1.2	0.6 - 1.8	1.2	0.6 - 1.8
BA (K/μL)	0.0	0.0 - 0.3		--Elev
NRBC (K/μL)				
NE (%)	66.5	58.5 - 74.5		
NE/MO/BA (%)			74.0	66.0 - 82.0
LY (%)	15.4	10.4 - 20.4	15.4	10.4 - 20.4
MO (%)	7.2	4.9 - 9.5		--Elev
EO (%)	10.6	5.6 - 15.6	10.6	5.6 - 15.6
BA (%)	0.3	0.0 - 2.3		--Elev
NRBC (%)				
RBC (M/μL)	6.27	5.98 - 6.56	6.27	5.98 - 6.56
HB (g/dL)	15.2	14.6 - 15.8	15.2	14.6 - 15.8
HCT (%)	37.0	32.9 - 41.3	37.0	32.9 - 41.3
MCV (fL)	59.1	55.1 - 63.1	59.1	55.1 - 63.1
MCH (pg)	24.3	22.3 - 26.5	24.3	22.3 - 26.5
MCHC (g/dL)	41.2	35.4 - 48.1	41.2	35.4 - 48.1
RDW (%)	21.5	19.5 - 23.5	21.5	19.5 - 23.5
PLT (K/μL)	400	335 - 465	400	335 - 465

**DREW SCIENTIFIC, INC. 14100 NW 57th Court Miami Lakes, FL 33014 • USA**

### IMPORTANT

Follow mixing instructions from reverse side:  
Wipe threads on bottle and cap before closing

**FOR INVITRO DIAGNOSTIC USE**

**PROTECT FROM FREEZING**

STORE AT 2°C - 8°C (35°F - 46°F)