



DREW SCIENTIFIC, INC.

MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: HORSE Product Number: 600063, 600071

Lot 2007 Expiration Date: 10/15/2017

For use with HEMAVET® and MASCOT™ Hematology Analyzers

Parameter	Hemavet 800,850,850FS,950,950LV 1500,1700,3700, Mascot 850		Hemavet 550, Mascot 550	
	Target	Range	Target	Range
WBC (K/μL)	17.6	16.1 - 19.1	17.6	16.1 - 19.1
NE (K/μL)	8.7	7.2 - 10.4		
NE/MO/BA (K/μL)			9.7	7.6 - 12.1
LY (K/μL)	6.9	5.5 - 8.5	6.9	5.5 - 8.5
MO (K/μL)	0.8	0.4 - 1.4		--Elev
EO (K/μL)	0.9	0.0 - 2.0	0.9	0.0 - 2.0
BA (K/μL)	0.2	0.0 - 1.0		--Elev
NRBC (K/μL)				
NE (%)	49.5	44.5 - 54.5		
NE/MO/BA (%)			55.35	47.4 - 63.4
LY (%)	39.3	34.3 - 44.3	39.3	34.3 - 44.3
MO (%)	4.8	2.5 - 7.1		--Elev
EO (%)	5.3	0.3 - 10.3	5.3	0.3 - 10.3
BA (%)	1.0	0.0 - 5.1		--Elev
NRBC (%)				
RBC (M/μL)	9.08	8.69 - 9.47	9.08	8.69 - 9.47
HB (g/dL)	16.8	16.2 - 17.4	16.8	16.2 - 17.4
HCT (%)	43.0	37.6 - 48.6	43.0	37.6 - 48.6
MCV (fL)	47.3	43.3 - 51.3	47.3	43.3 - 51.3
MCH (pg)	18.5	17.1 - 20.0	18.5	17.1 - 20.0
MCHC (g/dL)	39.1	33.4 - 46.3	39.1	33.4 - 46.3
RDW (%)	22.6	20.6 - 24.6	22.6	20.6 - 24.6
PLT (K/μL)	553	393 - 713	553	393 - 713

DREW SCIENTIFIC, INC. 14100 NW 57 Court Miami Lakes, FL 33014 • USA

IMPORTANT

Follow mixing instructions from reverse side:
Wipe threads on bottle and cap before closing

FOR INVITRO DIAGNOSTIC USE

PROTECT FROM FREEZING

STORE AT 2°C - 8°C (35°F - 46°F)