



# DREW SCIENTIFIC, INC.

## MULTI-TROL™

Multispecies Hematology Reference Control

Species Key: HORSE      Product Number: 600063, 600071

Lot 2047      Expiration Date: 10/15/2018

**For use with HEMAVET® and MASCOT™ Hematology Analyzers**

| Parameter       | Hemavet 800,850,850FS,950,950LV<br>1500,1700,3700, Mascot 850 |             | Hemavet 550, Mascot 550 |             |
|-----------------|---|-------------|-------------------------|-------------|
|                 | Target  | Range       | Target                  | Range       |
| WBC (K/μL)      | 18.5  | 17.0 - 20.0 | 18.5                    | 17.0 - 20.0 |
| NE (K/μL)       | 9.8   | 8.2 - 11.6  |                         |             |
| NE/MO/BA (K/μL) |   |             | 10.2                    | 8.0 - 12.6  |
| LY (K/μL)       | 7.3   | 5.9 - 8.9   | 7.3                     | 5.9 - 8.9   |
| MO (K/μL)       | 0.2   | 0.0 - 0.6   |                         | --Elev      |
| EO (K/μL)       | 1.0   | 0.1 - 2.1   | 1.0                     | 0.1 - 2.1   |
| BA (K/μL)       | 0.2   | 0.0 - 1.0   |                         | --Elev      |
| NRBC (K/μL)     |   |             |                         |             |
| NE (%)          | 53.1  | 48.1 - 58.1 |                         |             |
| NE/MO/BA (%)    |   |             | 55.03                   | 47.0 - 63.0 |
| LY (%)          | 39.6  | 34.6 - 44.6 | 39.6                    | 34.6 - 44.6 |
| MO (%)          | 0.9   | 0.0 - 3.2   |                         | --Elev      |
| EO (%)          | 5.3   | 0.3 - 10.3  | 5.3                     | 0.3 - 10.3  |
| BA (%)          | 1.0   | 0.0 - 5.1   |                         | --Elev      |
| NRBC (%)        |   |             |                         |             |
| RBC (M/μL)      | 8.62  | 8.23 - 9.01 | 8.62                    | 8.23 - 9.01 |
| HB (g/dL)       | 15.9  | 15.3 - 16.5 | 15.9                    | 15.3 - 16.5 |
| HCT (%)         | 41.7  | 36.5 - 47.2 | 41.7                    | 36.5 - 47.2 |
| MCV (fL)        | 48.4  | 44.4 - 52.4 | 48.4                    | 44.4 - 52.4 |
| MCH (pg)        | 18.4  | 16.9 - 20.0 | 18.4                    | 16.9 - 20.0 |
| MCHC (g/dL)     | 38.0  | 32.3 - 45.1 | 38.0                    | 32.3 - 45.1 |
| RDW (%)         | 23.2  | 21.2 - 25.2 | 23.2                    | 21.2 - 25.2 |
| PLT (K/μL)      | 551   | 391 - 711   | 551                     | 391 - 711   |

**DREW SCIENTIFIC, INC. 14100 NW 57 Court Miami Lakes, FL 33014 • USA**

### IMPORTANT

Follow mixing instructions from reverse side:  
Wipe threads on bottle and cap before closing

**FOR INVITRO DIAGNOSTIC USE**

**PROTECT FROM FREEZING**

STORE AT 2°C - 8°C (35°F - 46°F)